

Clinical Viral Transport Medium



MGIN Viral Collection Kit

Made in Korea

- Optimized swab for sample collection
- Convenient self-collection, storage and transport
- Easy to use for anyone with provided How-To-Use
- Simultaneous collection of nasopharyngeal and oral sample
- Compatible with various nucleic acid extraction kits



**DOH FDA RITM
VALIDATED**



**FDA U.S. FOOD & DRUG
ADMINISTRATION**

- *100 Tests per box
- *Room temperature storage



MAIN COMPOSITION	CAS No	Content (%)
Guanidine thiocyanate	593-84-0	35.0
Triton X-100	9002-93-1	3.9
Citrate sodium	68-04-2	1.0
8 - hydroxyquinoline	148-24-3	0.1
Deionized water	7732-18-5	≥80

pH 6-8

What is COVID-19 Sample Collection Kit?

This kit is called Clinical VTM (Viral Transport Medium). Clinical VTM is intended for the collection and transport of clinical specimen containing COVID-19 (Corona Virus) and others.

Why 2 Swabs?

For initial diagnostic testing for COVID-19, Korea CDC (Korea Center for Disease Control and Prevention) recommends collecting and testing upper respiratory (Nasopharyngeal(NP) swab and Oropharyngeal(OP) Swab).

Official Representative Philippines:

Medical@wwlgroup.biz / www.wwlifelink.com

(632) 5310-9184 loc.321/316 . 0917 135 4416 / 0977 820 8800